

Anti- Human B-cell receptor CD22/SIGLEC2



Product Name	Anti-Human B-cell receptor CD22/SIGLEC2
Catalog No	IT-300-125M1
Source	Mouse monoclonal antibody, clone 6G8C11
Isotype	IgG
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human B-cell receptor CD22/SIGLEC2 (IT-300-125p).
Applications	WB
Epitope	Linear
Specificity	Reacts with Human B-cell receptor CD22/SIGLEC2 protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

