

# Anti-Spike Protein of SARS-CoV-2 (CoVID-19/Wuhan), Rabbit Polyclonal



Product Name	Anti-Spike Protein of SARS-CoV-2 (CoVID-19/Wuhan)
Catalog No	IT-002-030
Source	Rabbit Polyclonal Antibody
Isotype	Rabbit IgG
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Plasmid vector expressing Spike protein of SARS-CoV-2 (CoVID-19/Wuhan).
Applications	WB, ELISA
Specificity	Reacts with Spike Protein S1 from Coronavirus SARS-CoV-2 (CoVID-19/Wuhan)(IT-002-030p), and with Spike Protein RBD domain from Coronavirus SARS-CoV-2 (CoVID-19/Wuhan)(IT-002-031p).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

