

Anti- Spike NTD (SARS-CoV-2/Wuhan)



Product Name	Anti- Spike NTD (SARS-CoV-2/Wuhan)
Catalog No	IT-002-039M1
Source	Mouse monoclonal antibody, clone 2F11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing S1 (SARS-CoV-2/Wuhan)(IT-002-030p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Bind to Spike NTD domain (SARS-CoV-2/Wuhan) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

