

## Anti-Spike RBD (CoVID-19/Wuhan)



Product Name	Anti-Spike RBD (CoVID-19/Wuhan)
Catalog No	IT-002-031M5
Source	Mouse Monoclonal Antibody, clone 1D11
Isotype	Mouse IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Spike RBD (CoVID-19/Wuhan)(IT-002-031p).
Applications	ELISA, FACS, etc.
Specificity	Reacts with CoV-19 Spike RBD, S1 and Spike (CoVID-19/Wuhan).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

