

Anti- Nucleocapsid protein (Borna disease virus /strain V)



Product Name	Anti- Nucleocapsid protein (Borna disease virus /strain V)
Catalog No	IT-028-002M3
Source	Mouse monoclonal antibody, clone 3D2F10
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Nucleocapsid protein (Borna disease virus /strain V) (IT-028-002p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Nucleocapsid protein (Borna disease virus /strain V) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

