

## Anti- Human CSF1R/Macrophage colony-stimulating factor 1 receptor



Product Name	Anti- Human CSF1R/Macrophage colony-stimulating factor 1 receptor
Catalog No	IT-300-123M1
Source	Mouse monoclonal antibody, clone 1A9H2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human CSF1R (IT-300-123p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Human CSF1R as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

