



p29 (GRA7)(Toxoplasma gondii)

Product Name	p29 (GRA7)(Toxoplasma gondii)
Catalog No	IT-011-053p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	p29 (GRA7)(Toxoplasma gondii) (aa 27-180; GenBank# O00933) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc.
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	

