



## Anti- Human RAET1L / ULBP6

Product Name	Anti-Human RAET1L / ULBP6
Catalog No	IT-300-149M5
Source	Mouse monoclonal antibody, clone 3G5G8
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human RAET1L / ULBP6 (IT-300-149p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Human RAET1L / ULBP6 as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

