


# GP (Ebola/Zaire)



Product Name	GP (Ebola/Zaire)
Catalog No	IT-014-003p-histag
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Ebola [Zaire] glycoprotein (aa 1-647) (Genbank No. AAN37507) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: GP1/GP2 (Ebola/Zaire) proteins

Note: