

# Anti-CP204L/p30 (ASFV/Pig/HLJ/2018)



---

Product Name	Anti-CP204L/p30 (ASFV/Pig/HLJ/2018)
Catalog No	IT-029-003M3
Source	Mouse monoclonal antibody, clone 5F6E6
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing CP204L/p30 (ASFV/Pig/HLJ/2018) (IT-029-003p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with CP204L/p30 (ASFV/Pig/HLJ/2018) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

