

# Anti-Glycoprotein G2 (Machupo mammarenavirus)



Product Name	Anti-Glycoprotein G2 (Machupo mammarenavirus)
Catalog No	IT-008-011M2
Source	Mouse monoclonal antibody, clone 1C9G5
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Glycoprotein G2 (Machupo mammarenavirus)(IT-008-011p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Glycoprotein G2 (Machupo mammarenavirus).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

