


HA(Δ TM)(A/Victoria/4897/2022)(H1N1)



Product Name	HA(Δ TM)(A/Victoria/4897/2022)(H1N1)
Catalog No	IT-003-00120 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	HA(Δ TM)(A/Victoria/4897/2022)(H1N1)(aa 18-529) with 6xHis tagged at C-terminus. Most in trimer.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	\geq 95%
Storage	Keep it at 4 $^{\circ}$ C if used within a month. For long term storage, split it into small aliquots and keep at -80 $^{\circ}$ C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein Seq:

DTLCIGYHANNSTDTVDTVLEKNVTVTHSVNLLEDKHNGLKCKLRGVAPHLGQCNIAGWILGNPECESLSTARSWSYIVETSNSDN
GTCYPGDFINYEELREQLSSVSSFERFEIFPKTSSWPNHSDNGVTAACSHAGARSFYKNLIWLKKGKSYPKINQTYINDKGKEVLVL
WGIHHPPTITDQESLYQNADAYVFGTISRYSKFKFPEIAARPKVRDQAGRMNYYWTLVEPGDKITFEATGNLVAPRYAFTMEKEAG
SGIISDTPVHDCNATCQTPEGAINSTLFPQNVHPITIGKCPKYVRSTKLRRLATGLRNVPSIQSRGLFGAIAAGFIEGGWTGMVDGWYGYH
HQNDQGGSGYAADLKSTQNAIDKITNKVNSVIEKMNTQFTAVGKEFNHLEKRIENLNKKVDDGFLDVWTYNAELLVLENERLTLDYH
DSNVKNLYEKVRHQLKNNAKEIGNGCFEFYHKCDNTCMESVKNGTYDYPKYSEEAKLNREKIDGVKLDSTRIYQ