

# Anti-Human Interleukin-33 / IL-33



Product Name	Anti-Human Interleukin-33 / IL-33
Catalog No	IT-300-160M3
Source	Mouse monoclonal antibody, clone 1D12G5
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human Interleukin-33 / IL-33 (IT-300-160Ep).
Applications	WB, ELISA, IP, IF, etc.
Epitope	Linear
Specificity	Reacts with Human Interleukin-33 / IL-33 as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

