

# Anti-EEV glycoprotein A35R [Monkeypox virus MPXV]



|                |  |
|----------------|--|
| Product Name   | Anti-EEV glycoprotein A35R [Monkeypox virus MPXV]  |
| Catalog No     | IT-012-015M1   |
| Source         | Mouse monoclonal antibody, clone 1G3B5   |
| Isotype        | IgG1   |
| Formulation    | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.  |
| Immunogen      | DNA vaccine expressing EEV glycoprotein A35R [Monkeypox virus MPXV] (IT-012-015p).   |
| Applications   | WB, ELISA, IP, IF etc.   |
| Epitope        | Linear   |
| Specificity    | Reacts with EEV glycoprotein A35R [Monkeypox virus MPXV] as tested.  |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

