

Anti-GP1 [LASSA Virus]



Product Name	Anti-GP1 [LASSA Virus]
Catalog No	IT-008-002M3
Source	Mouse monoclonal antibody, clone 1F12
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing GP1 [LASSA Virus] (IT-008-002p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with GP1 [LASSA Virus] as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

