

# Anti-A56(WR181)(Vaccinia Virus)



---

Product Name	Anti-A56(WR181)(Vaccinia Virus)
Catalog No	IT-012-006M5
Source	Mouse monoclonal antibody, clone 2D8E10
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing A56(WR181)(Vaccinia Virus) (IT-012-006p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with A56(WR181)(Vaccinia Virus) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

