

# Env (SRV-1)



Product Name	Env (SRV-1)
Catalog No	IT-019-006p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Envelope protein (Simian retrovirus 1/SRV-1) (a.a. 23-527) (GenBank Accession No. P04027) with 6xHis tag at C-terminal. Majority are trimers.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95%, based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

