


# VP1(Norovirus GII isolate Hu/GII.17[P17]/NewYork0024/2024/USA)



Product Name	VP1(Norovirus GII isolate Hu/GII.17[P17]/NewYork0024/2024/USA)
Catalog No	IT-013-005Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	6xHis N-terminal tagged VP1(Norovirus GII isolate Hu/GII.17[P17]/NewYork0024/2024/USA) (GenBank Accession No. XSF82115; a.a. 1-538).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: VP1 (Norovirus GII isolate  
Hu/GII.17[P17]/NewYork0024/2024/USA)